

# SAFETY MD-TRANSPORTER® USER GUIDE



# COPYRIGHT

Copyright 2021 Technimount E.M.S. Holding Inc. (Technimount), an affiliated company of Technologies CGC Inc. All Rights Reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language in any form or by any means without the written permission of Technimount or its affiliate companies.

#### **Trademark Attributions**

Safety MD-Transporter is a registered trademark of Technologies CGC Inc. Stryker, PowerPro-XT and Mx-Pro are registered trademarks from Stryker.

Technimount does not have any commercial relationship with Strkyer. Stryker does not promote or sell Technimount products.

In order to ensure continuous improvements to our products and services, we invite you to communicate your comments to our Customer Service Department at customerservice@technimount.com.

#### **OFFICE ADDRESS**

TECHNIMOUNT E.M.S. HOLDING INC. 3770 Jean Marchand Street, Suite 100-C Quebec (QC) G2C 1Y6 Canada

info@technimount.com www.technimount.com

T + 1 888.639.2758 F + 1 855.339.6351

#### NOTE

For any issues with your Technimount product, its components, or for any technical questions during the installation or maintenance, please send an email to <a href="mailto:techsupport@technimount.com">techsupport@technimount.com</a>.

Please have the serial number of your Technimount product available (as shown in the figure below) when calling Technimount Customer Service or Technical Support. Include the serial number in all written communications.

SERIAL NUMBER LOCATION




FIGURE 1. SAFETY MD-TRANSPORTER POST

# TABLE OF CONTENTS

COPYRIGHT	
CONTACT INFORMATION	3
INTRODUCTION	6
Purpose of the Document	6
Product Description	6
Intended use of the Product	6
Symbols and Definitions	
Warning / Caution / Note	
IDENTIFICATION OF COMPONENTS	
IV Pro Adapter	g
MICRO BASE	
Safety MD-Transporter	_
Specifications	10
NOTICE	
Product Illustration	-
SUMMARY OF SAFETY PRECAUTIONS	
Warning/Caution/Note	14
INSTALLATION GUIDE	16
Removing Packaging for Installation	16
Tools Required for Installation	16
Cot/Stretcher Orientation Diagram	17
Preparing the Cot/Stretcher for Installation	18
REMOVING THE IV POLE	18
Safety MD-Transporter Installation	
PREPARING THE IV PRO ADAPTER PLATE	
Installing the Micro Base on the IV Pro Adapter Plate Position on the Cot/Stretcher	
PATIENT RIGHT AND LEFT INSTALLATION	
INSTALLING THE IV PRO ADAPTER	
INSTALLING THE IV POLE AND HOLDER ONTO THE IV PRO ADAPTER	_
Installing the Safety MD-Transporter Post into the Micro Base	28
OPERATION GUIDE	30
General Operating Guidelines	30
Operating Bracket with the Safety MD-Transporter	
Installing the Bracket onto the Safety MD-Transporter	
ROTATING THE BRACKET ON THE SAFETY MD-TRANSPORTER	
REMOVING THE BRACKET FROM THE SAFETY MD-TRANSPORTER	
SAFETY MD-TRANSPORTER PART NUMBERS	
PART NUMBERS FOR STRYKER POWERPRO COT/STRETCHERS	
SAFETY MD-TRANSPORTER WITH THE VERTICAL MICRO BASE - VERSION 1	
MAINTENANCE GUIDE	
CLEANING THE SAFETY MD-TRANSPORTER	
CLEANING THE SAFETY MID-TRANSPORTER	
CLEANING T NOCESS	
REMOVAL OF IODINE COMPOUNDS	

# TABLE OF CONTENTS

Preventive Maintenance	
Inspection Process and Schedule	37
INSPECTION AND MAINTENANCE RECORD	
Training Record	39
Replacement Parts	40
WARRANTY	41
Warranty Policy	41
RETURN POLICY	41 42
Return Authorization	42
Damaged Merchandise	42
International Warranty Clause	42
Claim Process	43
ANNEX A	44
IV PRO ADAPTER- MX-PRO- SYSTEM OVERVIEW	
ASSEMBLY INSTRUCTIONS	•
SAE CERTIFICATION	51

# **SYMBOLS AND DEFINITIONS**

<u>SYMBOL</u>	DESCRIPTION
	Warning and Caution, special attention is required. Consult accompanying documents
	Safe working load symbol and Load balance symbol
	Pinch Point

# WARNING / CAUTION / NOTE

The word warning, caution, or note carry special meaning and should be carefully reviewed.

	WARNING Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.
	CAUTION Alerts the reader of a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device to avoid damage that may occur as a result of use or misuse.
NOTE	Provides special information about the product.

#### **IDENTIFICATION OF COMPONENTS**

#### **IV PRO ADAPTER**

- The IV Pro Adapter is a long black plate.
- The IV Pro Adapter is installed on the cot/ Stretcher side rail.
- It is needed to install the Micro Base, the Safety MD-Transporter and the IV Pole together on the cot/stretcher.



#### **MICRO BASE**

- The round base of the Safety MD-Transporter inserts into the Micro Base.
- The locking mechanism secures the Safety Transporter in place.
- The red quick release button is used to quickly release the device from the Micro Base.




#### SAFETY MD-TRANSPORTER

- The Safety MD-Transporter post is provided with a Micro Base.
- It has a round base at the bottom which inserts into the Micro Base.
- It comes with 2 blocks that are needed to install the IV pole. NOTE: the 2 blocks are different in size. The longer block must go on left side of the plate (opposite the pointed end).
- It can rotate 360° for easy access to the medical device at all times.
- The Micro Base on the top of the Safety MD-Transporter is used to hold the medical device.
- The Safety MD-Transporter can hold a maximum weight of 25 pounds.



# **PRODUCT ILLUSTRATION**



# INTRODUCTION

#### **PURPOSE OF THE DOCUMENT**

This User Guide is designed to assist you with the installation, operation and maintenance of the Technimount Safety MD-Transporter®.

#### PRODUCT DESCRIPTION

The Safety MD-Transporter is a mounting system that consists of a Micro Base, mounted vertically (version 1) on a post, installed on Technimount's IV Pro Adapter® on another Micro Base on the side frame of the cot/stretcher. It is designed to support and transport portable medical devices with a maximum weight of twenty-five (25) lbs. in pre-hospital and hospital environments. The Safety MD-Transporter can be removed and stored safely when not in use. When used in emergency vehicles with a cot/stretcher, the Safety MD-Transporter rotates 360° for easy access to the medical device during transport and loading/unloading of the patient on the cot/stretcher.

The Safety MD-Transporter is inserted into the Micro Base equipped with a locking mechanism and can accept any type of mounting bracket equipped with an Anti-Rotation Micro Disc on the back. The Safety MD-Transporter is available for the right or the left side of the cot/stretcher. The IV Pro Adapter is required to install the Safety MD-Transporter. Another version of the post (version 2) is also available with a horizontal Micro Base on the top to mount different equipment, according to the space available and other equipment installed around the mount.

#### INTENDED USE OF THE PRODUCT

The Safety MD-Transporter is a secure mounting system for portable medical devices, designed to be installed on Stryker PowerPro™ cot/stretchers. The mount is intended to support and transport medical equipment during a medical intervention of a traumatized, an ambulatory or non-ambulatory patient (including infants, children and adults).

Once the Safety MD-Transporter is installed, it facilitates the safe transport of smaller portable medical equipment (i.e. defibrillator and/or pumps) in emergency vehicles. The use of the Safety MD-Transporter was designed to be used on an ambulance cot/stretcher and is intended for front-line responders who have received the training necessary to operate the system in the field. It's designed to work with other products on a cot/stretcher, such as the Safety Arm System. The Safety MD-Transporter is intended for transport purposes only. It shouldn't carry more than one medical device. It is not intended for more than one unit or any additional equipment attached to the post.

The Safety MD-Transporter is compatible only with Technimount's mounting solutions that have an Anti-Rotation Micro Disc on the back for version 1 or a Micro Disc on the bottom for version 2.

# INTRODUCTION

# **TECHNICAL SPECIFICATIONS**

Product Name	Safety MD-Transporter	
Product Description	Bracket mounting system to secure medical devices on the left or right side of a stretcher during EMS and Critical Care Transport.	
Cot/Stretcher Compatibility*	Stryker® ProXT, PowerPro, model 6085-6086, 6500 and 6505	
Environment of Usage	EMS (Ground & Flight) – Hospital/Clinic, Military	
Product Section Material	<ul> <li>Safety MD-Transporter: Aluminum 6061-T6 Pole with Stainless Stee 304 Base</li> <li>Micro Base: Aluminum 6061-T6</li> <li>IV Pro Adapter: Aluminum 6061-T6</li> </ul>	
Product Finish	<ul> <li>Safety MD-Transporter: Silver-Anodized Finish</li> <li>Micro Base: Silver-Anodized Finish</li> <li>IV Pro Adapter: Black-Anodized Finish</li> </ul>	
Environmental Regulations	Complying with RoHS 3	
Temperature of Usage	- 35° C to 45° C	
Certification	Tested in Compliance with: SAE J3043 (26 G)	
Installation	Device Retained on Main Assembly with 2 Upper Arms with Red Quick Release Mechanisms and Lateral Blockers	
Attachment Options	Anti-Rotation Micro Disc, Bottom Micro Disc	
Transportation Compatibility System	Compatible with Technimount's:  • Bracket Pro Serie®  • Micro Base  • IV Pro Adapter™  • Safety Arm System™	
Features	<ul> <li>Secures medical devices directly on stretcher/cot during air and ground EMS and Critical Care Transport</li> <li>Installs on the left or right side of the stretcher</li> <li>Access to the medical device from all access (360° rotation)</li> <li>Red Quick Release Mechanism to easily insert and remove the bracket</li> <li>Resistant to biological fluids and Quaternary-type disinfectants</li> <li>Vertical and Horizontal Versions Available</li> </ul>	
Physical Overall Dimensions	Safety MD-Transporter (Vertical)  • Width 3.94 in, Depth 3.4 in, Height 15.96 in  Safety MD-Transporter (Horizontal)	
Websta	• Width 3.94 in, Depth 5.95 in, Height 11.13 in 6.8 lb	
Weight	<del></del>	
Maxiumum Weight Capacity	25 lb	
Part Number	<ul> <li>540-00-0L13 - Safety MD-Transporter (Vertical) with Micro Base</li> <li>540-00-0L14 - Safety MD-Transporter (Horizontal) with Micro Base</li> </ul>	
Others Models and Configuration Options	Contact Technimount EMS for more information or for other options at customerservice@technimount.com	

<sup>\*</sup> Product and medical equipment manufacturer names are trademarks™ or registered trademarks™ of their respective holders. Technimount does not have a commerical relationship with these medical equipment manufacturers.

# SUMMARY OF SAFFTY PRECAUTIONS

#### WARNING/CAUTION/NOTE

Carefully read and strictly follow all Warnings and Cautions listed in this manual. Servicing and maintenance should only be done by qualified personnel.



#### WARNING

- ✓ To avoid injury, hands should be clear of the Micro Base while installing the medical equipment.
- ✓ Improper usage of the Safety MD-Transporter can cause serious injury to the patient or operator. Operate only as described in this user guide.
- ✓ Do not modify the Safety MD-Transporter or any components of the system. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator.
- ✓ It is the responsibility of the operator of the cot/stretcher to ensure that the Safety MD-Transporter being used with Technimount products meets the installation specifications. Injury may result if non-compatible products are used with Technimount System products.
- ✓ Do not allow untrained staff to assist in the operation of the Safety MD-Transporter. Untrained technicians/staff can cause injury to the patient or themselves.
- ✓ Do not sit or ride on the pole, damage to the product could occur, and may result in injury to the patient or operator.
- ✓ Grasping the Safety MD-Transporter improperly can cause injury. Keep hands, fingers and feet away from moving parts.
- ✓ Operators must be able to safely lift the total weight of the Safety MD-Transporter (post and medical device), and any items attached to the medical device.
- ✓ Do not attach restraints to the post, the Micro Base, or any other place of the system. It may result in damage to the system or equipment and/or result in injury to the patient or operator.
- ✓ Do not place any additional equipment or objects on the Safety MD-Transporter, other than approved Technimount System device brackets for medical devices.
- ✓ Always avoid any improper pushing or violent manipulation, as it may result in an unbalanced cot/stretcher and may cause injury to the patient or operator and/or damage to the product.
- ✓ To avoid the risk of patient injury or equipment damage, ensure that you properly mount and secure the medical device to the Micro Base.
- ✓ To avoid the risk of patient injury or equipment damage, you must ensure that the disc attached to the bracket under the medical device is properly locked into the Micro Base base.
- ✓ Improper maintenance can cause injury or damage to the product. Maintain the Safety MD-Transporter as described in this user guide. Use only Technimount approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty.
- ✓ Failure to properly clean or dispose of corrosive products will increase the risk of premature damage and may cause injury to the patient or operator.
- Ensure proper hand placement on a high position on the post. Hands should be clear of Micro Base while loading and unloading the Safety MD-Transporter.
- ✓ Do not adjust or remove the Safety MD-Transporter while you are moving with the cot/stretcher.

# SUMMARY OF SAFFTY PRECAUTIONS

#### WARNING/CAUTION"NOTE (CONTINUED)



#### **CAUTION**

- ✓ Certified technical personnel familiar with the cot/stretcher's specifications should complete the installation of the Safety MD-Transporter. Refer to the cot/stretcher manufacturer's manual before installing the Safety MD-Transporter and ensure that the installation does not damage or interfere with the use of the cot/stretcher.
- ✓ Inspect regularly the Safety MD-Transporter and Micro Base for any issue and loose screws.

#### NOTE

- ✓ Some cleaning products are corrosive in nature and may cause damage to the product if used improperly. If such products are used to clean Technimount equipment, the equipment must be rinsed with clean water and thoroughly dried following the cleaning. Failure to properly rinse and dry the Safety MD-Transporter will leave a corrosive residue on the surface, possibly causing premature corrosion of critical components.
- ✓ Installation of the Safety MD-Transporter was designed to be placed on either side of the stretcher with the use of the **right** or the **left** Safety MD-Transporter (patient right or left). When ordering, it is important to identify which side the Safety MD-Transporter is to be located. The Safety MD-Transporter cannot be reversed as the black metal plate (IV Pro Adapter plate) has been designed for a specific side.
- A preventive maintenance program should be established for all Technimount equipment. Preventive maintenance may need to be performed more frequently, based on the level of use of the product.
- ✓ Failure to use authorized parts, and lubricants, etc., could cause damage to the Safety MD-Transporter and will void the warranty of the product.
- ✓ This user guide should be considered a permanent part of the Safety MD-Transporter and should remain with the product.
- ✓ Technimount continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your Safety MD-Transporter and this user guide.

If you have any questions, or for more information, please contact Technimount Customer Service at customerservice@technimount.com.

# Notice/Certification

The Safety MD-Transporter is tested in compliance with the essential requirements and other relevant provisions of the Directive 1999/5/EC and SAE J3043. For more information, please contact Technimount Customer Service (customerservice@technimount.com).

TECHNIMOUNT E.M.S. HOLDING INC. C/O: Regulatory Affairs 3770 Jean Marchand Street, Suite 100-C Quebec (QC) G2C 1Y6 Canada T+1888.639.2758 F+1855.339.6351 technimount.com info@technimount.com



#### **WARNING**

The end user and administrator are responsible to validate regulations and standards for safety in their region, in order to comply with applicable safety regulations. Technimount is not responsible to inform the end user and/or administrator of applicable legislation for safety in their area.

#### REMOVING PACKAGING FOR INSTALLATION

- 1. Unpack boxes and check all items for proper operation.
- 2. Ensure that all shipping and packaging materials have been properly removed from the product(s) prior to installation.
- 3. Identify all of the product components and hardware prior to starting the installation.
- 4. Some hardware or parts provided might not be applicable to the model of your cot/stretcher.
- 5. The Safety MD-Transporter must work properly before it is put into service.
- 6. Refer to the identification of components at the beginning of this user guide if needed.

#### **NOTE**

✓ Loose items or debris on the cot/stretcher can interfere with the operation of the Safety MD-Transporter and cot/stretcher. Keep the cot/stretcher and the floor clear.

#### TOOLS REQUIRED FOR INSTALLATION

Below are the basic tools required for installation. Other tools may also be needed. If you have any questions regarding the installation, or need assistance, please contact our Technical Support team at techsupport@technimount.com



FIGURE 2. ALLEN KEY 5/32 IN




FIGURE 3. PHILIPS PH2 SCREWDRIVER

#### NOTE

✓ Electrical tools are not recommended for screwing into the cot/stretcher frame, as there is a potential risk of damage to the threads.

# **COT/STRETCHER ORIENTATION DIAGRAM**



# Installation Guide

#### PREPARING THE COT/STRETCHER FOR INSTALLATION

#### REMOVING THE IV POLE

#### **NOTE**

- ✓ If the cot/stretcher has an IV pole located on the side you want to install the Safety MD-Transporter, you will need to remove the IV Pole.
- ✓ To keep the IV pole on the same side as the cot/stretcher, you will need to install it on the IV Pro Adapter.
- ✓ You will need to validate the location of the installation before removing the IV pole hardware. The IV Pro Adapter is designed for either the right or left side.
- ✓ Technimount has two (2) Safety MD-Transporter models, one with a vertical Micro Base (version 1) and the other with a horizontal Micro Base (version 2).
- ✓ The installation procedure provided in this user guide is for **Patient Right** installation. The procedure is the same for both right and left sides.
- ✓ You might not have to perform these steps if you do not have an IV pole installed.
- ✓ Please validate if your cot has the "knee Gatch" feature or not; It will dictate the exact IV Pro Adapter model, for the proper installation on the cot.
- ✓ Please refer to the cot/stretcher manufacturer's manual for more information about installation/re-installation of the IV pole.

# PREPARING THE COT/STRETCHER FOR INSTALLATION (CONTINUED)

#### REMOVING THE IV POLE (CONTINUED)

1. To access screws that hold the IV pole and the holder, lift the IV pole to its upper position.



- 2. Remove all four (4) screws.
- 3. Remove IV Pole and holder.
- 4. Keep the IV Pole and holder and set aside to install on the IV Pro Adapter plate in a later step.



#### SAFETY MD-TRANSPORTER INSTALLATION

#### **NOTE**


✓ For installation of the IV Pro Adapter on Styker Cot Mx-Pro, go to Annex A.

#### PREPARING THE IV PRO ADAPTER PLATE

1. Install the 2 blocks provided with Safety MD-Transporter onto the IV Pro Adapter plate.

#### **NOTE**


✓ Ensure to install the short block on the side of the IV Pro Adapter plate with pointed end.



- 2. Using Philips PH2 screwdriver, tighten all 4 screws (2 on each block).
- 3. Insert screws from the bottom of the IV Pro Adapter plate.

#### NOTE

- ✓ Installation of the blocks is performed underneath of the plate, in order to have the blocks on the upper side once installed.
- 4. IV Pro Adapter plate with the 2 blocks installed.






#### INSTALLING THE MICRO BASE ON THE IV PRO ADAPTER PLATE


1. Five (5) screws are provided to install the Micro Base onto the IV Pro Adapter plate.



2. Align with the five (5) screw holes.



- 3. Add a drop of Loctite #243 to the screws before installation.
- 4. Tighten screws, while alternating screwheads.



POSITION ON THE COT/STRETCHER
Patient Right and Left Installation



#### CAUTION

✓ It is the responsibility of the Safety MD-Transporter operator to ensure that the system being used on the cot/ stretcher meets the installation specifications. Injury may result if a non-compatible cot is used.

#### **NOTE**

- ✓ When ordering the Safety MD-Transporter, it is important to specify on which side the it will be installed on, as the System cannot be installed on the opposite side. There is a right and left Safety MD-Transporter designed for each side.
- ✓ The Safety MD-Transporter cannot be switched from one side to the other. You need to order a system according to the side that is needed.
- ✓ The location for the installation of the IV Pro Adapter is the same on patient right or left sides.




FIGURE 4. SAFETY MD-TRANSPORTER - INSTALLATION POSITION ON PATIENT RIGHT AND LEFT SIDES

#### **INSTALLING THE IV PRO ADAPTER**

#### **NOTE**

✓ During the installation process, Technimount System recommends the use of Loctite #242 to avoid premature loosening of screws during normal use.



The Loctite™ Brand shown above is provided as a suggestion only. Technimount and its subsidiaries are not related to, nor does it commercially endorse this brand.

#### **NOTE**

- ✓ Adhere to your organization's safety protocols and procedures for the following installation steps.
- 1. Identify the location on the cot/stretcher where you need to install the IV Pro Adapter.

#### **NOTE**

- ✓ Make sure the IV Pro Adapter will correspond to the feature of the stretcher/cot with or without the knee gatch. The knee gatch feature will require two (2) recesses on the IV Pro Adapter plate, the one without knee gatch will have only one (1) long recess to obtain the proper fit on the stretcher/cot.
- 2. Pull and hold the Safety Pin (red round handle) while putting a pair of pliers to hold the pin open.





#### INSTALLING THE IV PRO ADAPTER (CONTINUED)

3. When pliers are secured in place on the screw to hold the Safety Pin open, turn the red round handle counter clockwise to unscrew and remove the red round handle.



4. Remove the red round handle.



5. Locate all 4 screw holes on the stretcher frame.



#### **INSTALLING THE IV PRO ADAPTER (CONTINUED)**


- 6. Align IV Pro Adapter plate with the 4 screw holes on the cot/stretcher
- 7. Put a drop of Locktite #242 on each screw.
- 8. Tighten all four (4) screws while alternating screwheads



9. Push pin from the inside of the clamp block outward and place a pair of pliers on the screw to hold open while you replace the red round handle.




10. Put a drop of Locktite #242 onto the screw before replacing the red round handle. This will prevent premature loosening of the screws during normal use.



# INSTALLING THE IV PRO ADAPTER (CONTINUED)

- 11. Replace the red round handle while turning clockwise.
- 12. Screw in tightly.



# INSTALLING THE IV POLE AND HOLDER ONTO THE IV PRO ADAPTER

1. Locate the 2 blocks on either side of the IV Pro Adapter plate.



- 2. Install the IV Pro Holder onto the block located near the Clamp Block.
- 3. Align with the two (2) screw holes.
- 4. Tighten both screws.



#### INSTALLING THE IV POLE AND HOLDER ONTO THE IV PRO ADAPTER (CONTINUED)

- 5. Install the IV pole onto the IV Pro Adapter plate second block
- 6. Align with the two (2) scres holes.
- 7. Use the same screws provided with the IV Pole or similar screws provided, with same threads and length.
- 8. Tighten both screws.



9. When not in use, lower the IV pole and clip it into its holder.



# INSTALLING THE SAFETY MD-TRANSPORTER POST INTO THE MICRO BASE

1. Lift the IV Pole upright.



2. Slide the round base of the Safety MD-Transporter to secure it onto the Micro Base.




# INSTALLING THE SAFETY MD-TRANSPORTER POST INTO THE MICRO BASE (CONTINUED)

3. Safety MD-Transporter installed on the Micro Base.



4. If the IV Pole is not needed, you can fold and place the IV Pole into its holder. It fits perfectly beside the clamp block located on the cot/stetcher.



#### **GENERAL OPERATING GUIDELINES**

- ✓ Use the Safety MD-Transporter only as described in this User Guide.
- ✓ Read all the labels and instructions on the product before using the Safety MD-Transporter.
- ✓ Always inspect the Micro Base and the Safety MD-Transporter for any issue or damage before every use.
- ✓ Installing or removing the Safety MD-Transporter requires a minimum of one (1) trained operator.
- ✓ Do not remove the Safety MD-Transporter while moving with the cot/stretcher.
- ✓ Stay with the patient and control the Safety MD-Transporter at all times.





#### WARNING

✓ Never let a patient grab, manipulate or hold the Safety MD-Transporter, or the Micro Base on the cot/stretcher because it may cause injury to the operator or patient. Secure the Safety MD-Transporter as recommended while the patient is on the cot/stretcher.

# **OPERATING BRACKET WITH THE SAFETY MD-TRANSPORTER**

#### INSTALLING THE BRACKET ONTO THE SAFETY MD-TRANSPORTER

- Only Technimount brackets with an Anti-Rotation Micro Disc on the back can be inserted into the Micro Base located on the Safety MD-Transporter.
- 2. Ensure that your bracket has an Anti-Rotation Micro Disc on the back to install the medical device onto the Safety MD-Transporter.



3. While using both hands, slowly insert Anti-Rotation Micro Disc located on the back of the bracket into the Micro Base.



# OPERATING THE SAFETY MD-TRANSPORTER (CONTINUED)

#### INSTALLING THE BRACKET ONTO THE SAFETY MD-TRANSPORTER (CONTINUED)

4. A Technimount bracket with an Anti-Rotation Micro Disc on the back installed on the Safety MD-Transporter.



# ROTATING THE BRACKET ON THE SAFETY MD-TRANSPORTER


1. To rotate the Safety MD-Transporter with the medical device 360°, turn the Safety MD-Transporter clockwise or counter clockwise.



# **OPERATING THE SAFETY MD-TRANSPORTER (CONTINUED)**

#### REMOVING THE BRACKET FROM THE SAFETY MD-TRANSPORTER

- 1. To remove the bracket and medical device from the Safety MD-Transporter, press the red quick release button located on the back Micro Base.
- 2. While holding the red quick release button, lift the bracket upwards in a straight line.



# SAFETY MD-TRANSPORTER PART NUMBERS

#### PART NUMBERS FOR STRYKER POWERPRO COT/STRETCHERS

Some exceptions may apply for the PowerPro Stretcher, please confirm with Customer Service at customerservice@technimount.com prior to placing your order.

#### SAFETY MD-TRANSPORTER WITH THE VERTICAL MICRO BASE - VERSION 1

DESCRIPTION	<u>PART NUMBER</u>
The Safety MD-Transporter Post	540-10-OL13
2 Blocks for the IV Pro Adapter	543-00-UN
Micro Base with Quick Release Mechanism	120-20-UN
IV Pro Adapter – IV Pole Support – Patient Left Installation - Without Knee Gatch	511-00-UN-LFS - WHKG
IV Pro Adapter – IV Pole Support – Patient Left Installation - With Knee Gatch	511-00-UN-LFS - KG
IV Pro Adapter – IV Pole Support – Patient Right Installation - Without Knee Gatch	511-00-UN-RHS - WHKG
IV Pro Adapter – IV Pole Support – Patient Right Installation ¬With Knee Gatch	511-00-UN-RHS - KG

<sup>\*</sup> For available kits, please contact Customer Service at customerservice@technimount.com.

#### SAFETY MD-TRANSPORTER WITH THE HORIZONTAL MICRO BASE - VERSION 2

<u>DESCRIPTION</u>	PART NUMBER	
The Safety MD-Transporter Post (Angled Post)	540-10-OL14-ANG	
2 Blocks for the IV Pro Adapter	543-00-UN	
Micro Base with Quick Release Mechanism	120-20-UN	
IV Pro Adapter – IV Pole Support – Patient Left Installation - Without Knee Gatch	511-00-UN-LFS - WHKG	
IV Pro Adapter – IV Pole Support – Patient Left Installation - With Knee Gatch	511-00-UN-LFS - KG	
IV Pro Adapter – IV Pole Support – Patient Right Installation - Without Knee Gatch	511-00-UN-RHS - WHKG	
IV Pro Adapter – IV Pole Support – Patient Right Installation ¬With Knee Gatch	511-00-UN-RHS - KG	

<sup>\*</sup> For available kits, please contact Customer Service at customerservice@technimount.com.

#### **CLEANING THE SAFETY MD-TRANSPORTER**

#### **CLEANING PROCESS**

The Safety MD-Transporter can be cleaned by using a pressure washer. The unit may show some signs of oxidation or discoloration from continuous washing, however, degradation of the cot's performance or functionality will not be affected due to power washing if cleaning instructions are properly followed.

The material used to manufacture the Safety MD-Transporter is high-grade aluminum, it is very durable and resistant to corrosion, however high levels of sodium or a high concentration of sodium or acid may prematurely corrode the product. Always rinse with soft soap and water if exposed.

- ✓ Use of a power washer can accelerate removal of contaminants collected during the use of the product.
- ✓ Rinse with clean water.
- Towel dry all casters and interface points, inside the triangle head and clamp block.
- ✓ Allow all other components to air dry.
- Avoid over saturation and ensure that the product does not stay wet longer than the cleaner manufacturer's guidelines for proper disinfecting.



#### WARNING

✓ When cleaning, always use appropriate personal protection equipment (PPE) based on established protocols (e.g., gloves, eyewear, etc.).



#### **CAUTION**

- ✓ Do not steam clean or use ultrasonic cleaners on the system or any of its components.
- ✓ Do not immerse the metal parts/components in water.
- ✓ Maximum water temperature should not exceed 180°F/82°C
- ✓ Maximum water pressure should not exceed 1500 psi/103.5 BAR . If a pressure washer is being used to clean the unit, the pressure nozzle must be kept a minimum of 24 inches (61 cm) from the unit.

#### **CLEANING SOLUTIONS**

Phenolic type or quaternary type disinfectants can be used (excluding Virex® TB). Iodophor type disinfectants, Peracetic Acid, Chlorine, and especially Hydrogen Peroxide are not recommended for use because staining and corrosion may result.

Suggested cleaners for the Safety MD-Transporter surfaces and clamp box:

- ✓ Quaternary Cleaners (active ingredient ammonium chloride)
- ✓ Phenolic Cleaners (active ingredient o-phenyl phenol)
- ✓ Chlorinated Bleach Solution (3.25% less than 1-part bleach to 100 parts water)

#### **CLEANING THE SAFETY MD-TRANSPORTER (CONTINUED)**

#### **CLEANING SOLUTIONS (CONTINUED)**



#### WARNING

- ✓ Some cleaning products are corrosive in nature and may cause damage to the product if used improperly.
- ✓ If the products described above are used to clean the equipment, the Safety MD-Transporter must be rinsed with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the Safety MD-Transporter will leave a corrosive residue on the surface of the product, possibly causing premature corrosion of critical components.

#### REMOVAL OF IODINE COMPOUNDS

Use a solution of 1/4 Tablespoon of Sodium Thiosulfate in 1 pint / .5 litre of warm water to clean the stained area. Clean as soon as possible after staining occurs. If stains are not immediately removed, allow solution to soak or stand on the surface. Thoroughly rinse surfaces with clean water and allow to air dry before returning unit to service.



#### WARNING

✓ Failure to properly clean contaminated parts or components will increase the risk of exposure to blood borne pathogens and may cause injury to the patient or the operator.

#### PREVENTIVE MAINTENANCE

A preventive and regular maintenance program should be established for all Technimount equipment. Preventive maintenance may be required more frequently based on the usage level of the product. The Safety MD-Transporter requires regular preventive maintenance. Establish and follow a preventive maintenance schedule and keep records of preventive maintenance activities.



#### CAUTION

- ✓ Close attention should be given to safety features including, but not limited to both Micro Bases, Anti-Rotation Micro Disc, and post.
- ✓ Improper maintenance can cause injury or damage to the product. Maintain the product as described in this user guide.

#### LUBRICATION

The Safety MD-Transporter has been designed to operate without the need for lubrication. Due to an anodization of the metallic parts, using a silicone-base lubricant can be used but it is not necessary. If you're not sure about the product to use, please contact technical support or customer service.



#### **CAUTION**

✓ Failure to use authorized parts, lubricants, etc. could cause damage to the cot/stretcher and will void the warranty of the product.

# **INSPECTION PROCESS AND SCHEDULE**

#### PREVENTIVE MAINTENANCE PROGRAM

The following inspection routine and schedule is intended as a general guide for preventive maintenance of the Safety MD-Transporter. Factors such as weather, environment, geographical location, and individual usage will necessitate different maintenance. If you are unsure as to how to perform these maintenance inspections or at what interval to perform these inspections, please contact your Technical Support Team at techsupport@technimount.com. Check each routine and replace damaged or worn out parts if necessary.

Follow table below in accordance to the timetable, whichever comes first.

<u>ITEM</u>	ROUTINE	1 M OR 2 H	3 M OR 6 H	6 M OR 12 H	12 M OR 24 H
Micro Base	<ul> <li>✓ Inspect Mechanism</li> <li>✓ Inspect Red Quick Release Mechanism</li> <li>✓ Inspect General Functionalities</li> </ul>	X			
Anti- Rotation Back Disc	<ul> <li>✓ Inspect Screws that hold the Disc</li> <li>✓ Inspect the General Condition of the Disc</li> <li>✓ Regularly Clean the the Disc</li> </ul>	Х			
Safety MD- Transporter Post	<ul> <li>✓ Inspect Screws that hold the Micro Base</li> <li>✓ Inspect the General Condition of the Disc at the Base of the Post</li> <li>✓ Inspect the Assembly of the Post</li> <li>✓ Ensure There are no Loose Parts</li> </ul>	x			

# Maintenance Guide

# **INSPECTION AND MAINTENANCE RECORD**

PREVENTIVE MAINTENANCE PERFORMED	BY	<u>DATE</u>	<u>TIME</u>

Return To Table of Contents

# **TRAINING RECORD**

	TRAINING DATE		TRAINING
TRAINEE NAME	BASIC TRAINING	TRAINING UPDATE	TRAINING METHOD: USER GUIDE, IN-SERVICE, IN-CLASS, ETC.

38

# MAINTENANCE GUIDE

# **REPLACEMENT PARTS**

The parts and accessories listed are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this user guide may not be individually available for purchase. Please call Technimount Customer Service or Technical Support: +1.581-998-9820 or at customerservice@technimount.com or techsupport@technimount.com for availability and pricing.

DESCRIPTION	PART NUMBER
The Safety MD-Transporter Post (Version 1)	540-10-OL13
The Safety MD-Transporter Post (Angled Post) (Version 2)	540-10-OL14-ANG
2 Blocks for the IV Pro Adapter	543-00-UN
Micro Base with Quick Release Mechanism	120-20-UN
IV Pro Adapter – IV Pole Support – Patient Left Installation - Without Knee Gatch	511-00-UN-LFS - WHKG
IV Pro Adapter – IV Pole Support – Patient Left Installation - With Knee Gatch	511-00-UN-LFS - KG
IV Pro Adapter – IV Pole Support – Patient Right Installation - Without Knee Gatch	511-00-UN-RHS - WHKG
IV Pro Adapter – IV Pole Support – Patient Right Installation ¬With Knee Gatch	511-00-UN-RHS - KG
Screws for the Micro Base (Phillips)	544-00-UN
Screws for the IV Pro Adapter	543-00-UN



# **WARRANTY POLICY**

This statement constitutes TECHNIMOUNT's entire warranty policy with regards to the TECHNIMOUNT Products.

TECHNIMOUNT MAKES NO OTHER WARRANTY OR REPRESENTATION NEITHER EXPRESSED NOR IMPLIED, EXCEPT AS SET FORTH HEREIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL TECHNIMOUNT BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

TECHNIMOUNT E.M.S. Holding Inc. (TECHNIMOUNT) guarantees to the original "Purchaser" of the "Hardware" with which this "Limited Warranty" is included that the Hardware will be free from "Defects" in workmanship and materials under normal use for a "Warranty Period", one (1) year from the date that the Hardware was first purchased by the Purchaser. During the Warranty Period, the Hardware will be repaired or replaced according to the "Limited Warranty" without charge to the Purchaser for either parts or labor. The parts may be repaired or replaced with either new or refurbished parts. The product may be replaced for a new or refurbished product. For this Limited Warranty "refurbished" means parts and products, which have been returned to factory specifications. If the Hardware is repaired or replaced within the Warranty Period, the greater of the remaining warranty period will apply or three (3) months from the date of repair/replacement. If the Hardware is repaired or replaced after the Warranty Period has expired, the Warranty Period for the repair or replacement will expire three (3) months after the date of repair or replacement.

#### LIMITED RESPONSIBILITY AND WARRANTY

The Limited Warranty does not apply to normal wear that should result from normal use. It does not apply when the Hardware or any component is opened, disassembled, or repaired by someone not authorized by TECHNIMOUNT and does not cover repair or replacement of any Hardware or part thereof damaged by neglect, misuse, moisture, liquids, exposure to heat, accidents, abuse, and non- compliance with the instructions for installation and use provided with the Hardware. The Limited Warranty does not cover physical damage to the surface of the Hardware, and the decision to repair, replace or refuse the coverage is final and its sole discretion of TECHNIMOUNT, and is made without any compensation or obligation from TECHNIMOUNT.

The hardware, or product defined as a mounting or bracket system for clipping and attaching medical equipment is specifically designed to fill this requirement. Any other use will void the warranty and TECHNIMOUNT will not be held liable on any claim if the product is modified or adapted for use.

Technimount products are intended to retain a medical device in place in the case of a single crash impact. A Technimount product must NOT be used again if it was involved in a crash. The product MUST be replaced. If the Purchaser uses a Technimount product following a crash, it is at the Purchaser's own risk and Technimount will not be held liable.

#### INTERNATIONAL WARRANTY CLAUSE

This warranty reflects Canadian domestic policy. Warranty outside Canada may vary by country. Please contact TECHNIMOUNT Customer Service for more information.



#### **RETURN POLICY**

Mounting systems and bracket systems for portable medical devices may be returned up to 60 days of receipt if:

- The product is not what was originally ordered
- The product does not meet specifications with TECHNIMOUNT's technical sheets according to the customers
- The product is not compatible with the equipment on which it was to be installed

For any manufacturing defect, refer to the conditions of the Warranty Policy.

#### PRIOR TO 30 DAYS

- Product must be undamaged and in its original packaging
- The product return request must be provided in writing and it must be approved by TECHNIMOUNT prior to returning the product
- Returns will not be approved on a modified or damaged item
- Charges may apply if damaged or missing item
- Customer is responsible for a 10% restocking fee

## PRIOR TO 45 DAYS

- Product must be unused, undamaged and in its original packaging
- The product return request must be provided in writing and it must be approved by TECHNIMOUNT SYSTEM prior to returning the product
- Returns will not be approved on modified or damaged item
- Charges may apply if damaged or missing items
- Customer is responsible for a 25% restocking fee

#### PRIOR TO 60 DAYS

- Product must be unused, undamaged and in its original packaging
- The product return request must be provided in writing and it must be approved by TECHNIMOUNT prior to returning the product
- Returns will not be approved on a modified or damaged item
- Charges may apply if damaged or missing items
- Customer is responsible for a 30% restocking fee

## **RETURN AUTHORIZATION**

TECHNIMOUNT's customer service department is responsible for approving any merchandise return and will provide a Return Merchandise Authorization (RMA) number to be printed on any returned merchandise. TECHNIMOUNT reserves the right to charge shipping and restocking fees on returned items. SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.



#### DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receiving merchandise. DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT. Upon prompt notification, TECHNIMOUNT will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by TECHNIMOUNT within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Note: Claims for any short or broken shipment must be made within thirty (30) days of invoice.

## **CLAIM PROCESS**

Please follow the claim procedure below to return a defective product:
Fill in the RMA form online at https://www.technimount.com/service-support/service-support-overview/

#### REQUIRED INFORMATION

- Serial number
- Invoice number
- Identification of the issue
- Upload the pictures of the issue
- Submit the form

A confirmation of the receipt of the claim will be sent to you within 2 business days. If you have not received a confirmation Email by then, please contact techsupport@technimount.com.

If the RMA request is approved, TECHNIMOUNT will provide you with an RMA number for your claim. We will also include instructions for the replacement of your product. No product return will be accepted without the RMA number labeled on the outside of the box being returned.


#### FINDINGS AND CONCLUSION

- If following an investigation, we find that our product is defective, TECHNIMOUNT will take corrective actions and close the claim
- If following an investigation, we find that our product is NOT defective (misuse or abuse of the product) the product will not be covered by the warranty. Details of our findings and conclusions will be provided to you.

## **QUESTIONS ABOUT OUR POLICY**

If you have any questions regarding our return policy, please contact our Customer Service department at customerservice@technimount.com.

# IV PRO ADAPTER - MX-PRO - SYSTEM OVERVIEW



# IV PRO ADAPTER - MX-PRO - SYSTEM OVERVIEW (CONTINUED)

# <u>IDENTIFICATION OF THE COMPONENTS OF THE IV PRO ADAPTER</u>

<u>ITEM</u>	PART NUMBER	DESCRIPTION	QUANTITY
1	511-00-MXP-LFS	IV Pro Adapter Plate	1
2	530-10-MXP3-RHS	Safety Arm System Clamp Block - Mx-Pro	1
3	0004-163-000	Button Head Cap Screw (from Stryker Mx-Pro Assembly)	2
4	TCM-039-002	I.V. Clip Spacer	1
5	TCM-039-003	I.V. Pole Spacer	1
6	6080-110-035	I.V. Clip (from Stryker Mx- Pro Assembly)	1
7		I.V. Pole (from Stryker Mx- Pro Assembly)	1
8	TCM-999-012	Hex Drive BO Button Head Cap Screw 1/4-20 x 1"	2
9	TCM-VR-999-008	Hex Drive SS Button Head Cap Screw 1/4-20 x 7/16"	4
10	TCM-VF-999-185	Phillips SS Flat Head Screw 10-32 x 5/8"	4
11	120-20-UN	Micro Base	1
12	TCM-VF-999-173	Phillips SS Flat Head Screw 1/4-20 x 3/4" with Nylon Patch	5

# ANNEX A

## **ASSEMBLY INSTRUCTIONS**


# WARNING/CAUTION/NOTE




## **CAUTION**

- ✓ Only certified technical personnel familiar with the bracket installation and use should manipulate it. Always refer to Technimount's user guide for all instructions for installing and using any Technimount System mounting system.
- ✓ Inspect regularly all components on the equipment for any issue and loose screws, bolts and nuts.
- ✓ Never install the mount into or on other similar mounting systems or brackets without written confirmation by Technimount System as these different systems may be from other brands or models and may not be compatible with Technimount System.


1. Locate the Patient Left Outer Rail Extrusion.



2. Unscrew the I.V. Pole from the stretcher. Store the Stryker screws for future use.



3. Unscrew the I.V. Clip from the Outer Rail Extrusion. Keep the Stryker screws for Step #8.




4. Install the I.V. Clamp and the I.V. Pole Spacer on the IV Pro Adapter Plate using four (4) Phillips Stainless Steel Flat Head Screw 10-32 x 5/8".



## **CAUTION**

✓ Apply Loctite #243 on every screw.




5. Install the I.V. Clamp and the I.V. Pole in the Spacers using four (4) Hex Drive Stainless Steel Button Head Cap Screw 1/4-20 x 7/16".



## **CAUTION**

✓ Apply Loctite #243 on every screw.




6. Assemble the Safety Arm System Clamp Block on the Stretcher Outer Rail Extrusion according to the Safety Arm System user guide procedure.




## **CAUTION**

✓ Apply Loctite #243 on every screw.



7. Place the IV Pro Adapter on the Stretcher Outer Rail Extrusion. Align the screw holes of the plate with the threaded holes of the Stretcher Outer Rail Extrusion.




 Install the IV Pro Adapter to the Outer Rail Extrusion using two (2) screws removed at Step #3 and two (2) Hex Drive Black Oxide Button Head Cap Screw 1/4-20 x 1".




# CAUTION


✓ Apply Loctite #243 on every screw.



9. Install the Micro Base on the IV Pro Adapter using five (5) Phillips Stainless Steel Flat Head Screw 1/4-20 x 3/4" with Nylon Patch.



10. The system is ready.







# SAFETY AND FLEXIBILITY WHERE IT MATTERS MOST

# TECHNIMOUNT E.M.S. HOLDING INC.

3770 Jean Marchand Street, Suite 100-C Quebec (QC) G2C 1Y6 Canada T + 1 888.639.2758

F + 1 855.339.6351 technimount.com

info@technimount.com

