



**TECHNIMOUNT**  
EMS™

Technimount E.M.S. Holding Inc.

3770 Jean-Marchand Street, Suite 100-C,  
Québec, QC, G2C 1Y6,  
Canada

Phone: +1-581-700-6613  
Toll Free: 888-639-2758

[info@technimount.com](mailto:info@technimount.com) [www.technimount.com](http://www.technimount.com)

## Spacer Pro™ For Bracket Pro Serie™ mounts

### Features

- Lightweight
- Anodized coating
- Compact size
- Does not block the 360° rotation of swivel base
- Same holds
- Safe & Easy
- Immediate access
- Matches perfectly with our mounting systems

### Application of usage

- Crash Cart
- Emergency Work Space
- Stretcher
- Biomed Dept.

### Advantages

- Minimum space required for installation
- No additional holes needed
- Latex free
- Biological fluid and sodium resistant
- Rust resistant
- 1-year full warranty

Any medical equipment that has been inadequately installed in ambulances, other emergency medical service (EMS) vehicles or crash carts, can potentially pose a serious risk or become dangerous in the event that the vehicle makes frequent or sudden stops, turns or rushed departures. This situation can result in damaged expensive medical equipment or even inflict injuries on a patient, ambulance personnel or clinic staff.



Our option mount systems **Spacer Pro** are designed for a fast solution to raise or extend the surface base installation, with no interferences to the surface base. The features of this extension support are its very small size (weight), perfect fit to the standard base, with aligned holds for less time for installation. They can be installed in the same position or area to secure the mount and are offered at six (6) different heights or can be customized. The mounting extensions offer different heights for universal or standard surface base.

They have many characteristics; more specifically allow rotating up to 360 degrees with the base, fitting the shape and the holes for quick installation, without any compromise of use the surface base. The surface system extension will keep the features of the base and resist to an impact of many times the G- Force, in an accidental drop situation. To spare the expenses of costly repair or replacement of the equipment, and medical care and possible legal redress for injury, consider outfitting your crash cart, emergency workspace or biomed workshop with a deployment system for monitor and/or defibrillator.

### About Us

At Technimount EMS, we understand that EMS and Critical Care Transport teams need solutions that they can rely on. We design and develop quality mounting systems which are safe and flexible, allowing you to concentrate on what matters most - patient care.

For use with [Technimount System Bracket Pro Series™](#) mounting systems

## TECHNICAL SPECIFICATIONS

<b>Product name:</b>	<b>Spacer Pro™</b>	
<b>Compatibility model:</b>	Bracket Pro Serie - Standard or Universal model only	
<b>Environment of usage:</b>	EMS - Hospital - Clinic - Biomed Dept. - Military hospital	
<b>Material:</b>	Base: Aluminum high density with anodization process	
<b>Temperature of usage:</b>	- 35° C @ 45° C (Direct conduction)	
<b>Impact Resistance<sup>1</sup>:</b>	<b>Drop test 1.5 M and Same as Bracket Pro model</b>	
<b>No. of installation points:</b>	4 points under surface base aligned with the surface base	
<b>Dimensions (both parts assy - overall):</b>	7-1/2" With. X 8" Depth. X 3/4" <b>Height depending.</b> (in)	
<b>Rotation:</b>	Allows 360° degrees <b>with standard base</b>	
<b>Features:</b>	Compatible with <b>Standard and Extended Surface Base</b>	
<b>Weight:</b>	Total: 1.54 lb - 2.58 lb (depending on height)	
<b>Options:</b>	Spacer (height extension mount 3/8", 1/2", 3/4", 1", 1-1/4")	
<b>Part number<sup>2</sup>:</b>	825-00-SPC038	Spacer Pro - 3/8" height
<b>Part number</b>	825-00-SPC012	Spacer Pro - 1/2" height
<b>Part number</b>	825-00-SPC034	Spacer Pro - 3/4" height
<b>Part number:</b>	825-00-SPC100	Spacer Pro - 1" height
<b>Part number</b>	825-00-SPC114	Spacer Pro - 1-1/4" height
<b>Part number:</b>	825-00-SPC114-AN	Spacer Pro - 1-1/4" Inclined height

Counter Lip Height (in)	Surface Base Height (in)	Spacer Height Required (in)	Part No.
1/4"	1 1/4"	3/8"	825-00-SPC038
3/8"	1 1/4"	3/8"	825-00-SPC038
1/2"	1 1/4"	1/2"	825-00-SPC012
3/4"	1 1/4"	3/4"	825-00-SPC034
1"	1 1/4"	1"	825-00-SPC100
1 1/4"	1 1/4"	1 1/4"	825-00-SPC114
Inclined 1 1/4"	1 1/4"	Inclined 1 1/4"	825-00-SPC114-AN



### *Suggested installation*

Spacer Pro, Technimount EMS are registered trademarks of Technologies CGC Inc.  
 © 2021 Technimount E.M.S. Holding Inc. - All Rights Reserved - Printed in Canada. .