

## Arm Mount Pro adapter™

### Bracket mounting system for Arm installation

#### Characteristics

- Lightweight
- Compact size
- Quick Release and re-attachment
- Safe and Easy
- Immediate access
- Match with **Bracket Pro** series mount system from us.

#### Application of use

- ER Wall installation
- Arm extension
- Intensive care
- Arm support as GCX mounts and similar

#### Advantages

- Minimum space for installation
- For all kind of device
- Latex free
- Resist to sodium and the biological fluid
- Insensitive to rust
- Full warranty for 1 year

Any medical equipment that has been inadequately installed in ambulances, helicopters and other emergency medical service vehicles (EMS) and crash carts, can potentially pose a serious risk or become dangerous in the event that the vehicle makes frequent or sudden stops, turns or rushed departures. This situation can result in damaged expensive medical equipment or even inflict injuries on a patient, ambulance personnel or clinic staff.



Our mount systems “**Arm Mount Pro Adapter™**” are designed for specific installation on a wall using a flexible arm in a rail, for fast access to the unit, in times of emergency uses. The features of this support are the very small size (weight) and the specific positioning of

installation on a wall in emergency room or clinics. This new adapter is designed to be installed on a surface base using an arm support on a wall to give many positions to have access to the medical equipment. This brand-new system is designed to get the use of the medical equipment easier during a patient intervention. The Arm Mount Pro Adapter accepts the surface base **Bracket Pro** series (*Universal or short base only*). Perfect option for retrofit or to be installed on a GCX system using arm support. The base support is locked on the platform and the rotation feature is still active with this, however limited by the distance from the wall to the back of the unit. The system has been tested to meet the manufacturer requirements and provide with all the safety and security in hospital environment. The installation of the product is very simple and takes only just a few minutes, on the base.

The new **Arm Mount Pro Adapter** allows you to install the **Bracket Pro™** on wall, (with flexible and adjustable arm) instead of a working surface or crash cart. The **Arm Mount Pro Adapter™** comes with a standard (or extended) surface base and it can be adapted with all the Bracket Pro. The deployment of the unit is almost instantaneously and the system locked the unit after use.

**Technimount System Inc** sells, develops and commercialize deployment system for medical equipment that can be fit with different brand names of defibrillator, in hospital, EMS, or military environment.

For use with model [Technimount](#) Series « [Bracket Pro™](#) »

## **TECHNICAL SPECIFICATIONS**

<b>Product name:</b>	<b>Arm Mount Pro Adapter</b>
<b>Compatibility model:</b>	Bracket Pro series 20, 25, 30, 35, 40, 45 and 50
<b>Environment of use:</b>	EMS / Hospital – Clinic – Biomed shop – Military hospital -
<b>Material:</b>	Base: Aluminum high density w/ anodize process
<b>Temperature of use:</b>	- 35° C @ 45° C (Direct conduction)
<b>Impact Resistance<sup>1</sup>:</b>	<b>(Drop test 1.5 m)</b> According to the Bracket Pro model
<b>Nbr of installation point:</b>	6 points on surface base: 4 points on plate adapters
<b>Dimensions support:</b>	7-1/2" Width X 11" Depth X 5/16" Thick. (inch.)
<b>Dimensions plate adaptor:</b>	4-5/8" length X 4-13/16" width X 1/4" Thickness (inch.)
<b>Rotation:</b>	360° Surface base on the arm support <b>with locking system</b>
<b>Mechanical Installation:</b>	Sliding insert with locking pin and set screws
<b>Features:</b>	<u>Compatible with:</u> Standard or Extended surface base from <b>Technimount GCX mounting solutions</b> arms and supports
<b>Weight:</b>	Support Assy: 1.55 lbs. – 703 g Support Assy Standard version: 1.39 lbs – 589 g
<b>Options:</b>	Standard or Extended Surface base
<b>Part number<sup>2</sup>:</b>	601-00-WMA-GCX Adapter GCX arm with Extended base 602-00-WMA-GCX Adapter GCX arm with Standard base



*Installation suggested*

<sup>1</sup> *Technimount System Inc* works to meet and exceed all legal requirements in the deployment system for EMS vehicles, and hospital regulations according to the Provincial & Federal Government regulations, it may not meet your local requirements, please contact your administrator.

<sup>2</sup> Part number must be confirmed prior sending purchase order. *Technimount System inc* reserves the right to make changes on the product, documents and technical specifications without notice and responsibilities. DOCUMENTS PRINTED IN CANADA  
CGC Technologies – Division of Medical equipment & accessories – Registered brand 2013.