Technimount System

Administration
445 Papin Street, Office 103
Quebec City (Quebec) Canada
G1P 3T8

Telephone: +1-581-998-9820 Fax: 1-855-339-6351

Info@technimount.com

www.technimount.com

Bracket Pro Serie™ 30

Mounting System for Monitor / Defibrillator Physio-Control LP-12

Features

- > Lightweight
- Compact size
- Swivel base 360° with locking system
- Quick Release & Reattachment
- Safe & Easy
- Immediate access
- Matches with our other mounting systems

Application of usage

- Crash Cart
- Ambulance
- EMS Vehicles
- Stretcher/Cot
- Military Vehicles

Advantages

- Minimum space required for installation
- No vertical release
- Latex free
- Biological fluid resistant
- Rust resistant
- 1-year full warranty

Any medical equipment that has been inadequately installed in ambulances, other emergency medical service (EMS) vehicles or crash carts, can potentially pose a serious risk or become dangerous in the event that the vehicle makes frequent or sudden stops, turns or rushed departures. This situation can result in damaged expensive medical equipment or even inflict injuries on a patient, ambulance personnel or clinic staff.



Our locking systems *Bracket Pro SerieTM*30 are designed for fast access to the unit, with little effort in times of emergency uses. Features include lightweight design, incorporates a Quick Release & Reattachment mechanism and anodized finish of 15 microns for scratch and corrosion resistance. When inserted into the *Standard Surface Base*, it allows a

360-degree rotation of the medical device while securely locked.

To spare the expenses of costly repair or replacement of the equipment, medical care and possible legal redress for injury, consider outfitting your EMS vehicle or crash cart with a deployment system for monitor and/or defibrillator. The mounting systems offer different configurations for most external defibrillators on the market. The deployment of the unit is nearly instantaneous when engaging the release button on the base. The complete assembly has been tested for G-Force tolerances and compliant with SAE standards.

About Us

Technimount System creates mounts and brackets for portable medical devices such as defibrillators, ventilators, and pumps, that adapt to different brands of medical equipment for OEM ambulance manufacturers, EMS users, and hospitals.

For use with Physio-Control Defibrillator LIFEPAK® 12

TECHNICAL SPECIFICATIONS

Product name:	Bracket Pro Serie™ 30
Compatibility model:	Physio-Control - LIFEPAK® 12 Monitor / Defibrillator
Environment of usage:	EMS / Hospital - Clinic - BioMed Dept Military hospital
Material:	Base: High density aluminum w/ anodization process Device section: High density aluminum w/ anodization process
Temperature of usage:	- 35° C @ 45° C (Direct conduction)
Impact Resistance1:	Compliant with industry requirements
No. of installation points:	4 points on surface base: 9 points on device mount
Dimensions (both parts assy- overall):	12-1/8" Wid. X 8" Depth. X 3-1/4" Thickness overall (in)
Rotation:	360° degrees with locking system (Standard Surface Base) 15° degrees both sides with locking system (Extended Surface Base)
Features:	 With or without Carry case protector Compatible with Technimount Stretcher Safety Arm System™ Compatible with Extended Surface Base Free access to all connectors and accessories (top)
Weight:	Total (2 parts) 2.9 lb - device section: 1.30 lb
Options:	Extended Surface Base (longer) Spacer Pro (height extension mount 3/8", ½", ¾", 1", 1 ¼") Converter Pro (surface base cover for work space)
Locking mechanism resistance:	Tested 20,000 + "press & release" activations without any reduction in use
Part number ² :	200-00-PC12 Bracket Pro Serie™ 30 Complete Assembly
Part number:	200-10-PC12 Bracket Pro Serie™ 30 Bracket Device Only
Part number:	100-20-UN Standard Surface Base Only





Suggested installation

¹ *Technimount System* works to meet and exceed all legal requirements in the deployment system for EMS vehicles, according to the Provincial & Federal Government regulations. It may not meet your local requirements, please contact your administrator.

² Part number must be confirmed prior tp sending purchase order. Technimount System reserves the right to make changes to the product, documents and technical specifications without notice and responsibilities. DOCUMENTS PRINTED IN CANADA.
CGC Technologies – Division of Medical equipment & accessories – Registered brand 2013.
PHYSIO-CONTROL® and LIFEPAK12® are trademarks and registered trademarks of PHYSIO-CONTROL.